


Hindfoot Fusion Nail



FracTech | Trauma System


Division of **Absolute** Medical Devices Pvt Ltd


INDICATION

- Severe deformity of the foot or ankle
- Arthritis
- Instability and skeletal defects following tumor resection
(including, but not limited to, neuro-osteoarthropathy like Charcot's foot)
- Avascular necrosis of the talus
- Failed joint replacement or ankle fusion
- Nonunions of Distal tibial fracture
- Osteoarthritis
- Rheumatoid arthritis and pseudoarthrosis


CONTRAINDICATIONS

- Dysvascular limbs
- Active infections
- Insufficient plantar pad





Hindfoot Fusion Nail



Ø 10mm (RIGHT)

Code	Length(mm)
FT32.10.150R	150mm
FT32.10.240R	240mm
FT32.10.180R	180mm

(LEFT)

Code	Length(mm)
FT32.10.150L	150mm
FT32.10.240L	240mm
FT32.10.180L	180mm

Ø 12mm (RIGHT)

Code	Length (mm)
FT32.12.150R	150mm
FT32.12.240R	240mm
FT32.12.180R	180mm

(LEFT)

Code	Length(mm)
FT32.12.150L	150mm
FT32.12.240L	240mm
FT32.12.180L	180mm

Ø 13mm (RIGHT)


Code	Length(mm)
FT32.13.150R	150mm
FT32.13.240R	240mm
FT32.13.180R	180mm

(LEFT)

Code	Length(mm)
FT32.13.150L	150mm
FT32.13.240L	240mm
FT32.13.180L	180mm

Ø 5.0 mm INTERLOCKING BOLT

with T25 Star Drive Recess



Code	Length (mm)
FT50LB.26	26mm
FT50LB.28	28mm
FT50LB.30	30mm
FT50LB.32	32mm
FT50LB.34	34mm
FT50LB.36	36mm
FT50LB.38	38mm
FT50LB.40	40mm
FT50LB.42	42mm


Code	Length (mm)
FT50LB.44	44mm
FT50LB.46	46mm
FT50LB.48	48mm
FT50LB.50	50mm
FT50LB.52	52mm
FT50LB.54	54mm
FT50LB.56	56mm
FT50LB.58	58mm
FT50LB.60	60mm

Code	Length (mm)
FT50LB.65	65mm
FT50LB.70	70mm
FT50LB.75	75mm
FT50LB.80	80mm
FT50LB.85	85mm
FT50LB.90	90mm
FT50LB.95	95mm
FT50LB.100	100mm

Ø 6.0 mm INTERLOCKING BOLT

with T25 Star Drive Recess

Code	Length (mm)	Code	Length (mm)	Code	Length (mm)
FT60LB.26	26mm	FT60LB.44	44mm	FT60LB.65	65mm
FT60LB.28	28mm	FT60LB.46	46mm	FT60LB.70	70mm
FT60LB.30	30mm	FT60LB.48	48mm	FT60LB.75	85mm
FT60LB.32	32mm	FT60LB.50	50mm	FT60LB.80	80mm
FT60LB.34	34mm	FT60LB.52	52mm	FT60LB.85	85mm
FT60LB.36	36mm	FT60LB.54	54mm	FT60LB.90	90mm
FT60LB.38	38mm	FT60LB.56	56mm	FT60LB.95	95mm
FT60LB.40	40mm	FT60LB.58	58mm	FT60LB.100	100mm
FT60LB.42	42mm	FT60LB.60	60mm		



SPIRAL BLADE



Code	Length (mm)	Code	Length (mm)
FT61.SB.40	40mm	FT61.SB.75	75mm
FT61.SB.45	45mm	FT61.SB.80	80mm
FT61.SB.50	50mm	FT61.SB.85	85mm
FT61.SB.55	55mm	FT61.SB.90	90mm
FT61.SB.60	60mm	FT61.SB.95	95mm
FT61.SB.65	65mm	FT61.SB.100	100mm
FT61.SB.70	70mm		

END CAP, FOR HINDFOOT FUSION NAIL






Code
FT62.EC.20

END CAP




Code
FT63.EC.30

Hindfoot Fusion Nail Implant Box



Hindfoot Fusion Nail Instrument Set





FracTech | Trauma System

Division of **Absolute** Medical Devices Pvt Ltd



ABSOLUTE MEDICAL DEVICES (OPC) PVT LTD

E 39 1st Floor Part-2 SMA Industrial Area, GT Karnal Road, Delhi-110033

absolutemedicaldevices@gmail.com

+91 9810518801