


# Advance Proximal Femur Nailing System





# About us

Fractures can often be unexpected and difficult to avoid.  
When faced with a fracture, seeking the appropriate treatment is crucial  
For proper healing and Recovery..

As an Orthopedic implant manufacturer, our company plays a  
Vital role in the medical industry by designing, manufacturing,  
And distributing orthopedic implants used in the treatment of various  
Musculature conditions and injuries.  
These implants are crucial in restoring mobility, improving quality of life,  
And aiding in the recovery pf patients.

## Who We ARE?

We are an orthopedic implant manufacturer and have been supplying  
Implants to doctors and hospitals in the North India region.  
As an orthopedic implant manufacturer, our role is crucial in providing  
Healthcare professionals with the necessary tools to perform  
Orthopedic surgeries and improve the lives of patients.

By supplying orthopedic implants, we contribute to restoring mobility,  
Relieving pain, and enhancing the overall quality of life for patients  
Who require orthopedic procedures.  
Our products are used in various orthopedic surgeries, such as trauma  
Repairs, helping patients regain their mobility and functionality.




## Advance Proximal Femur Nail Long

This nail is intended for fixation of Low and extended subtrochanteric fractures, Ipsilateral trochanteric fractures, Combination of fractures (trochanteric area/shaft) and Pathological fractures.

The Blade helps in compacting the cancellous bone which provides a better hold especially in case of any osteoporotic bone.


### Proximal Locking Options



Ø10.5mm Helical Blade


Ø10.5mm Lag Screw

### Long Nail Distal Locking Option



Static

Dynamic



Ø4.9mm Locking Bolt

## Advance Proximal Femur Nail Short

This nail is intended for fixation of petrochanteric fractures, intertrochanteric fracture and High subtrochanteric fractures.

### Proximal Locking Options



Ø10.5mm Helical Blade



Ø10.5mm Lag Screw

### Long Nail Distal Locking Option




Static Oblique



Static (90°)



Dynamic



# Advance Proximal Femur Nail Short Titanium

Size	Item Code
170mm	FT24.PF09.17
200mm	FT24.PF09.20
240 mm	FT24.PF09.24
170mm	FT24.PF10.17
200mm	FT24.PF10.20
240 mm	FT24.PF10.24
170mm	FT24.PF11.17
200mm	FT24.PF11.20
240 mm	FT24.PF11.24
170mm	FT24.PF12.17
200mm	FT24.PF12.20
240 mm	FT24.PF12.24



## Advance Proximal Femur Nail Long Titanium

Left	Right
<b>FT24.PF09.34L</b>	<b>FT24.PF09.34R</b>
<b>FT24.PF09.36L</b>	<b>FT24.PF09.36R</b>
<b>FT24.PF09.38L</b>	<b>FT24.PF09.38R</b>
<b>FT24.PF09.40L</b>	<b>FT24.PF09.40R</b>
<b>FT24.PF09.42L</b>	<b>FT24.PF09.42R</b>
<b>FT24.PF09.44L</b>	<b>FT24.PF09.44R</b>

Left	Right
<b>FT24.PF10.34L</b>	<b>FT24.PF10.34R</b>
<b>FT24.PF10.36L</b>	<b>FT24.PF10.36R</b>
<b>FT24.PF10.38L</b>	<b>FT24.PF10.38R</b>
<b>FT24.PF10.40L</b>	<b>FT24.PF10.40R</b>
<b>FT24.PF10.42L</b>	<b>FT24.PF10.42R</b>
<b>FT24.PF10.44L</b>	<b>FT24.PF10.44R</b>

Left	Right
<b>FT24.PF11.34L</b>	<b>FT24.PF11.34R</b>
<b>FT24.PF11.36L</b>	<b>FT24.PF11.36R</b>
<b>FT24.PF11.38L</b>	<b>FT24.PF11.38R</b>
<b>FT24.PF11.40L</b>	<b>FT24.PF11.40R</b>
<b>FT24.PF11.42L</b>	<b>FT24.PF11.42R</b>
<b>FT24.PF11.44L</b>	<b>FT24.PF11.44R</b>

Left	Right
<b>FT24.PF12.34L</b>	<b>FT24.PF12.34R</b>
<b>FT24.PF12.36L</b>	<b>FT24.PF12.36R</b>
<b>FT24.PF12.38L</b>	<b>FT24.PF12.38R</b>
<b>FT24.PF1Ø2.40L</b>	<b>FT24.PF1Ø2.40R</b>
<b>FT24.PF1Ø2.42L</b>	<b>FT24.PF1Ø2.42R</b>
<b>FT24.PF1Ø2.44L</b>	<b>FT24.PF1Ø2.44R</b>



## Advance Blade for Proximal Femur Nail Titanium

Size	Item Code
70mm	FT25.BL70
75mm	FT25.BL75
80mm	FT25.BL80
85mm	FT25.BL85
90mm	FT25.BL90
95mm	FT25.BL95
100mm	FT25.BL100
105mm	FT25.BL105
110mm	FT25.BL110
115mm	FT25.BL115
120mm	FT25.BL120



## Advance Kephalikos Screw for Proximal Femur Nail Titanium


Size	Item Code
70mm	FT26.CE70
75mm	FT26.CE75
80mm	FT26.CE80
85mm	FT26.CE85
90mm	FT26.CE90
95mm	FT26.CE95
100mm	FT26.CE100
105mm	FT26.CE105
110mm	FT26.CE110





## Ø4.9mm Locking Bolts Self Tapping Titanium

Size	Item Code
26mm	FT49LB.26
28mm	FT49LB.28
30mm	FT49LB.30
32mm	FT49LB.32
34mm	FT49LB.34
36mm	FT49LB.36
38mm	FT49LB.38
40mm	FT49LB.40
42mm	FT49LB.42
44mm	FT49LB.44
46mm	FT49LB.46
48mm	FT49LB.48
50mm	FT49LB.50
52mm	FT49LB.52
54mm	FT49LB.54
56mm	FT49LB.56
58mm	FT49LB.58
60mm	FT49LB.60
65mm	FT49LB.65
70mm	FT49LB.70
75mm	FT49LB.75
80mm	FT49LB.80



## End Cap Advance Proximal Femur Nail Titanium

Size	Item Code
0mm	FT27.PF.EC00
5mm	FT27.PF.EC15
10mm	FT27.PF.EC10
15mm	FT27.PF.EC15



## Inner Screw for Advance Kephalikos Screw Titanium

Size	Item Code
0mm	FT27.PF.IS5



## End Cap for Advance Kephalikos Screw Titanium

Size	Item Code
0mm	FT27.PF.CS5





**Distributed & Marketed by:**

**Absolute** Health Solutions

Email : [ah.solutionsindia@gmail.com](mailto:ah.solutionsindia@gmail.com)

Ph. 9871299988,9319191251

