

Master Tibia Nailing System



FracTechTM | Trauma System
Division of **Absolute** Medical Devices Pvt Ltd



About us

Fractures can often be unexpected and difficult to avoid.
When faced with a fracture, seeking the appropriate treatment is crucial
For proper healing and Recovery..

As an Orthopedic implant manufacturer, our company plays a
Vital role in the medical industry by designing, manufacturing,
And distributing orthopedic implants used in the treatment of various
Musculoskeletal conditions and injuries.

These implants are crucial in restoring mobility, improving quality of life,
And aiding in the recovery pf patients.

Who We ARE?

We are an orthopedic implant manufacturer and have been supplying
Implants to doctors and hospitals in the North India region.

As an orthopedic implant manufacturer, our role is crucial in providing
Healthcare professionals with the necessary tools to perform
Orthopedic surgeries and improve the lives of patients.

By supplying orthopedic implants, we contribute to restoring mobility,
Relieving pain, and enhancing the overall quality of life for patients
Who require orthopedic procedures.

Our products are used in various orthopedic surgeries, such as trauma
Repairs, helping patients regain their mobility and functionality.

Indications

The Master Tibia Nail is Indicated for fractures in the tibial shaft as well as for metaphyseal and Certain intra-articular fractures off the Tibial Head and the pilon Tibial.


Master Tibia Nailing System offers Multiple Locking Options Which allows use of one Nailing System for many Tibia Bone Fracture patterns .

Proximal Locking




4.9mm Locking Bolts
in all Proximal holes

Distal Locking



3.9mm Locking Bolts in Ø8mm Nail
4.9mm Locking Bolts in dia 9mm, 10mm, 11mm



- Anterior - Posterior Locking Option
- Medial to Lateral Locking Option

Master Tibia Interlocking Nail Titanium

Code	Dia x Length
FT779.8.280	Ø8.0mm x 280mm
FT779.8.300	Ø8.0mm x 300mm
FT779.8.320	Ø8.0mm x 320mm
FT779.8.340	Ø8.0mm x 340mm
FT779.8.360	Ø8.0mm x 360mm
FT779.8.380	Ø8.0mm x 380mm

Code	Dia x Length
FT779.9.280	Ø9.0mm x 280mm
FT779.9.300	Ø9.0mm x 300mm
FT779.9.320	Ø9.0mm x 320mm
FT779.9.340	Ø9.0mm x 340mm
FT779.9.360	Ø9.0mm x 360mm
FT779.9.380	Ø9.0mm x 380mm


Code	Dia x Length
FT779.10.280	Ø10.0mm x 280mm
FT779.10.300	Ø10.0mm x 300mm
FT779.10.320	Ø10.0mm x 320mm
FT779.10.340	Ø10.0mm x 340mm
FT779.10.360	Ø10.0mm x 360mm
FT779.10.380	Ø10.0mm x 380mm

Code	Dia x Length
FT779.11.280	Ø11.0mm x 280mm
FT779.11.300	Ø11.0mm x 300mm
FT779.11.320	Ø11.0mm x 320mm
FT779.11.340	Ø11.0mm x 340mm
FT779.11.360	Ø11.0mm x 360mm
FT779.11.380	Ø11.0mm x 380mm




Ø4.9mm Locking Bolts Self Tapping Titanium

Size	Item Code
24mm	FT49LB.24
26mm	FT49LB.26
28mm	FT49LB.28
30mm	FT49LB.30
32mm	FT49LB.32
34mm	FT49LB.34
36mm	FT49LB.36
38mm	FT49LB.38
40mm	FT49LB.40
42mm	FT49LB.42
44mm	FT49LB.44
46mm	FT49LB.46
48mm	FT49LB.48
50mm	FT49LB.50
52mm	FT49LB.52
54mm	FT49LB.54
56mm	FT49LB.56
58mm	FT49LB.58
60mm	FT49LB.60
65mm	FT49LB.65
70mm	FT49LB.70
75mm	FT49LB.75
80mm	FT49LB.80



Ø3.9mm Locking Bolts Self Tapping Titanium

Size	Item Code
24mm	FT39LB.24
26mm	FT39LB.26
28mm	FT39LB.28
30mm	FT39LB.30
32mm	FT39LB.32
34mm	FT39LB.34
38mm	FT39LB.38
40mm	FT39LB.40
42mm	FT39LB.42
44mm	FT39LB.44
46mm	FT39LB.46
50mm	FT39LB.48
50mm	FT39LB.50
52mm	FT39LB.52
54mm	FT39LB.54
56mm	FT39LB.56
58mm	FT39LB.58
60mm	FT39LB.60





End Cap for Master Tibia Nail Titanium

Size	Item Code
0mm	FT30.EC.00
5mm	FT30.EC.05
10mm	FT30.EC.10
15mm	FT30.EC.15



P 206.00 Instrument Set for Master Tibia Nailing System





Distributed & Marketed by:

Absolute Health Solutions

Email : ah.solutionsindia@gmail.com

Ph. 9871299988,9319191251