

Multi Locking Proximal Humerus Nailing System



About us

Fractures can often be unexpected and difficult to avoid.
When faced with a fracture, seeking the appropriate treatment is crucial
For proper healing and Recovery..

As an Orthopedic implant manufacturer, our company plays a
Vital role in the medical industry by designing, manufacturing,
And distributing orthopedic implants used in the treatment of various
Musculoskeletal conditions and injuries.
These implants are crucial in restoring mobility, improving quality of life,
And aiding in the recovery pf patients.

Who We ARE?

We are an orthopedic implant manufacturer and have been supplying
Implants to doctors and hospitals in the North India region.
As an orthopedic implant manufacturer, our role is crucial in providing
Healthcare professionals with the necessary tools to perform
Orthopedic surgeries and improve the lives of patients.

By supplying orthopedic implants, we contribute to restoring mobility,
Relieving pain, and enhancing the overall quality of life for patients
Who require orthopedic procedures.
Our products are used in various orthopedic surgeries, such as trauma
Repairs, helping patients regain their mobility and functionality.

Indications

Multi Locking Proximal Humerus

Nail (Short)

This Multilocking Proximal Humerus Nail is indicated for fractures of the Proximal Humerus, including:


- 2 -Part Surgical neck Fractures
- 3 -Part Fractures
- 4 -Part Fractures

Multi Locking Proximal Humerus

Nail (Long)

This Multilocking Proximal Humerus Nail is indicated for fractures of the Proximal Humerus, including:

- 2 -Part Surgical neck Fractures
- 3 -Part Fractures
- 4 -Part Fractures
- Fractures of the Proximal Humerus with diaphyseal Extension
- Diaphyseal Fractures of the Humeral shaft
- Impending Pathological Fractures



Features

This Nailing System includes several innovative features designed specifically for the clinical needs of humeral nailing

- Straight Nail with central insertion point
- Unique screw-in-screw option, for improved Stability where Needed
- Ascending screw option to provide support to the calcar
- Multiplanar distal locking for less Implant Toggling
- Fully Targeted Distal Locking for the Short Nail

Multi Locking Proximal Humerus Nail

Size	Left	Right
Ø7.0mm x 160mm	FT774.7.160L	FT774.7.160R
Ø7.0mm x 200mm	FT774.7.200L	FT774.7.200R
Ø7.0mm x 220mm	FT774.7.220L	FT774.7.220R
Ø7.0mm x 240mm	FT774.7.240L	FT774.7.240R
Ø7.0mm x 260mm	FT774.7.260L	FT774.7.260R
Ø7.0mm x 280mm	FT774.7.280L	FT774.7.280R
Ø7.0mm x 300mm	FT774.7.300L	FT774.7.300R



Size	Left	Right
Ø8.0mm x 160mm	FT775.8.160L	FT775.8.160R
Ø8.0mm x 200mm	FT775.8.200L	FT775.8.200R
Ø8.0mm x 220mm	FT775.8.220L	FT775.8.220R
Ø8.0mm x 240mm	FT775.8.240L	FT775.8.240R
Ø8.0mm x 260mm	FT775.8.260L	FT775.8.260R
Ø8.0mm x 280mm	FT775.8.280L	FT775.8.280R
Ø8.0mm x 300mm	FT775.8.300L	FT775.8.300R

Size	Left	Right
Ø8.5mm x 160mm	FT776.85.160L	FT776.85.160R
Ø8.5mm x 200mm	FT776.85.200L	FT776.85.200R
Ø8.5mm x 220mm	FT776.85.220L	FT776.85.220R
Ø8.5mm x 240mm	FT776.85.240L	FT776.85.240R
Ø8.5mm x 260mm	FT776.85.260L	FT776.85.260R
Ø8.5mm x 280mm	FT776.85.280L	FT776.85.280R
Ø8.5mm x 300mm	FT776.85.300L	FT776.85.300R

4.5mm FracTech™ Multi Locking Screw for Humerus Nail

Size	Item Code
20mm	FT777.ML20
22mm	FT777.ML22
24mm	FT777.ML24
26mm	FT777.ML26
28mm	FT777.ML28
30mm	FT777.ML30
32mm	FT777.ML32
34mm	FT777.ML34
36mm	FT777.ML36
38mm	FT777.ML38
40mm	FT777.ML40
42mm	FT777.ML42
44mm	FT777.ML44
46mm	FT777.ML46
48mm	FT777.ML48
50mm	FT777.ML50
52mm	FT777.ML52
54mm	FT777.ML54
56mm	FT777.ML56
58mm	FT777.ML58
60mm	FT777.ML60



Ø3.9mm Locking Bolts

Size	Item Code
24mm	FT39LB.24
26mm	FT39LB.26
28mm	FT39LB.28
30mm	FT39LB.30
32mm	FT39LB.32
34mm	FT39LB.34
36mm	FT39LB.36
38mm	FT39LB.38
40mm	FT39LB.40
42mm	FT39LB.42
44mm	FT39LB.44
46mm	FT39LB.46
48mm	FT39LB.48
50mm	FT39LB.50
52mm	FT39LB.52
54mm	FT39LB.54
56mm	FT39LB.56
58mm	FT39LB.58
60mm	FT39LB.60



3.5mm FracTech™ Locking Head Screw



Size	Item Code
10mm	FT331.LH10
12mm	FT331.LH12
14mm	FT331.LH14
16mm	FT331.LH16
18mm	FT331.LH18
20mm	FT331.LH20
22mm	FT331.LH22
24mm	FT331.LH24
26mm	FT331.LH26
28mm	FT331.LH28
30mm	FT331.LH30
32mm	FT331.LH32
34mm	FT331.LH34
36mm	FT331.LH36
38mm	FT331.LH38
40mm	FT331.LH40
45mm	FT331.LH45
50mm	FT331.LH50
55mm	FT331.LH55
60mm	FT331.LH60
65mm	FT331.LH65
70mm	FT331.LH70
75mm	FT331.LH75
80mm	FT331.LH80
85mm	FT331.LH85
90mm	FT331.LH90



End Cap for Humerus Nail Titanium

Size	Item Code
0mm	FT90.221EC
2mm	FT90.222EC
5mm	FT90.225EC
10mm	FT90.210EC
15mm	FT90.215EC



P 209.00 Instrument Set for Multi Locking Humerus Nailing System



Distributed & Marketed by:

Absolute Health Solutions

Email : ah.solutionsindia@gmail.com

Ph. 9871299988, 9319191251